

# SEM EBIC Analysis Software: Technical Report

Panteha Farhadi

September 26, 2025



## Abstract

This report documents a comprehensive Python software package for analyzing Scanning Electron Microscopy (SEM) Electron Beam Induced Current (EBIC) measurements. The software processes multi-frame TIFF images with embedded metadata, extracts current maps, detects p-n junctions, and calculates diffusion lengths and depletion widths using exponential fitting algorithms. The graphical interface enables interactive visualization and analysis.

# Contents

<b>1</b>	<b>Introduction to EBIC Physics</b>	<b>4</b>
1.1	Voltage Dependence of Depletion Width . . . . .	4
<b>2</b>	<b>Software Architecture Overview</b>	<b>5</b>
<b>3</b>	<b>Core Components</b>	<b>5</b>
3.1	Metadata Extraction ( <b>Metadata</b> Class) . . . . .	5
3.1.1	Spatial Calibration Parameters . . . . .	5
3.1.2	Amplifier and Signal Processing Parameters . . . . .	5
3.1.3	Beam Parameters . . . . .	6
3.1.4	Scan Parameters . . . . .	6
3.1.5	Detector Parameters . . . . .	6
3.1.6	Image Acquisition Parameters . . . . .	6
3.1.7	Sample and Stage Parameters . . . . .	6
3.1.8	Instrument Identification . . . . .	6
3.1.9	Calibration Parameters . . . . .	6
3.2	Current Map Conversion ( <b>CurrentMap</b> Class) . . . . .	7
3.2.1	Pixel to Voltage Conversion . . . . .	7
3.2.2	Voltage to Current Conversion . . . . .	7
3.2.3	SEMProcessor Class . . . . .	8
3.3	Junction Detection Algorithms . . . . .	9
3.3.1	JunctionAnalyzer Class Overview . . . . .	9
3.3.2	Preprocessing with Bilateral Filtering . . . . .	9
3.3.3	Canny Edge Detection with Otsu Thresholding . . . . .	9
3.3.4	Post-Processing for Junction Modeling . . . . .	10
3.3.5	Coordinate Transformation . . . . .	11
3.3.6	Validation Metrics . . . . .	11
3.4	Diffusion Length Extraction Algorithms . . . . .	12
3.4.1	DiffusionLengthExtractor Class Overview . . . . .	12
3.4.2	Iterative Fitting Architecture . . . . .	12
3.4.3	Stage 1: Head Search for Stable Starting Point . . . . .	12
3.4.4	Stage 2: Tail Truncation for Noise Mitigation . . . . .	13
3.4.5	Iterative Convergence Loop . . . . .	14
3.4.6	Physical Model and Parameter Extraction . . . . .	14
3.4.7	Parameter Initialization for Fitting . . . . .	15
3.5	Interactive GUI for Visualization and Analysis . . . . .	16
3.5.1	Interactive GUI ( <b>SEMViewer</b> Class) . . . . .	16
3.6	Iteration over Multiple Images for Bias Analysis . . . . .	16
3.6.1	Overview . . . . .	16
3.6.2	Global Image Alignment . . . . .	16
3.6.3	Line Alignment . . . . .	17
3.6.4	Final Data Extraction and Validation . . . . .	18
<b>4</b>	<b>Analysis Workflow</b>	<b>20</b>
4.1	Step 1: Data Loading and Conversion . . . . .	20
4.2	Step 2: Junction Detection . . . . .	20
4.3	Step 3: Profile Analysis . . . . .	20

4.4	Step 4: Exponential Fitting . . . . .	20
<b>5</b>	<b>Software Workflow Overview</b>	<b>21</b>
5.1	GUI-Based Software Workflow . . . . .	21
5.1.1	TIFF File Processing . . . . .	21
5.1.2	Main Analysis Interface . . . . .	23
5.1.3	Image Visualization Tools . . . . .	23
5.1.4	Junction Analysis Workflow . . . . .	24
5.2	Parameter Extraction and Results . . . . .	27
5.2.1	Iterative Stabilization Algorithm for Diffusion Length Extraction	27
5.2.2	Algorithm Workflow . . . . .	27
5.2.3	Convergence Criteria Enhancement . . . . .	28
5.2.4	Enhancement of Physical Model: Bi-exponential Diffusion . . . .	28
5.2.5	Significance . . . . .	29
5.3	Pixel matching . . . . .	29
5.3.1	Profile Matching Precision . . . . .	31
<b>6</b>	<b>Future Development Plans</b>	<b>32</b>
6.1	Enhanced Edge Detection Based on Current Map . . . . .	32
6.2	Software Version Integration . . . . .	32
6.3	Convergence Criteria Enhancement . . . . .	32
6.3.1	Enhancement of Physical Model: Bi-exponential Diffusion . . . .	32
6.4	Advanced Pixel Matching Algorithms . . . . .	32
6.5	Pixel Matching and Profile Matching Precision . . . . .	32

# 1 Introduction to EBIC Physics

Electron Beam Induced Current (EBIC) is a characterization technique used in semiconductor physics to study minority carrier diffusion lengths and junction properties. When an electron beam strikes a semiconductor sample, it generates electron-hole pairs. In the presence of electric fields (such as near p-n junctions), these carriers separate and generate a measurable current. The current profile across a junction follows an exponential decay, from which key parameters can be extracted:

- **Diffusion length (L):** Characteristic distance minority carriers travel before recombination
- **Depletion width (W):** Region around the junction where electric fields dominate
- $1/\lambda$ : Inverse decay constant obtained from exponential fitting of EBIC profiles

The relationship is described by:

$$I(x) = I_0 \exp\left(-\frac{|x - x_j|}{L}\right) \quad (1)$$

where  $x_j$  is the junction position and  $L$  is the diffusion length.

## 1.1 Voltage Dependence of Depletion Width

The depletion width  $W$  in a p-n junction depends on the applied bias voltage  $V$  through the following relationship:

$$W(V) = \sqrt{\frac{2\epsilon}{q} \left( \frac{1}{N_A} + \frac{1}{N_D} \right) (V_{bi} - V)} \quad (2)$$

where:

- $W(V)$  is the voltage-dependent depletion width
- $\epsilon$  is the permittivity of the semiconductor material
- $q$  is the elementary charge ( $1.602 \times 10^{-19}$  C)
- $N_A$  and  $N_D$  are the acceptor and donor concentrations, respectively
- $V_{bi}$  is the built-in potential of the junction
- $V$  is the applied bias voltage (negative for reverse bias)

For a one-sided abrupt junction where  $N_A \gg N_D$  or vice versa, the equation simplifies to:

$$W(V) = \sqrt{\frac{2\epsilon}{qN} (V_{bi} - V)} \quad (3)$$

where  $N$  is the doping concentration on the lightly doped side.

## 2 Software Architecture Overview

The software consists of four main components:

1. TIFF metadata extraction and current conversion
2. Junction detection algorithms
3. Diffusion length extraction via exponential fitting
4. Interactive GUI for visualization and analysis
5. Iteration over multiple images for bias analysis using pixel matching to track profiles accurately

## 3 Core Components

### 3.1 Metadata Extraction (Metadata Class)

```
1 class Metadata:
2     def __init__(self, xml_path):
3         self.xml_path = xml_path
4         self.data = self._parse_metadata()
5
6     def _parse_metadata(self):
7         ns = {'rdf': 'http://www.w3.org/1999/02/22-rdf-syntax-ns#',
8              'image': 'http://ns.pointelectronic.com/Image/1.0/',
9              'efa': 'http://ns.pointelectronic.com/EFA/1.0/'}
10        # ... XML parsing logic
```

Listing 1: Metadata extraction from XML

This class extracts the following comprehensive set of SEM parameters from XML metadata:

#### 3.1.1 Spatial Calibration Parameters

- PixelSizeX, PixelSizeY: Physical size of each pixel (meters)
- PixelUnit: Unit of measurement for pixel dimensions

#### 3.1.2 Amplifier and Signal Processing Parameters

- Contrast, Brightness: Image contrast and brightness settings
- EffectiveAmpGain, Amplification: Amplifier gain factors
- InputOffset, OutputOffset: Signal offset corrections
- InverseEnabled: Whether signal inversion is applied

### 3.1.3 Beam Parameters

- **BeamVoltage:** Electron beam acceleration voltage (kV)
- **BeamCurrent:** Electron beam current (nA)
- **WorkingDistance:** Distance from objective lens to sample (mm)

### 3.1.4 Scan Parameters

- **ScanRotation:** Image rotation angle (degrees)
- **ScanSpeed:** Beam scanning speed
- **ScanDwellTime:** Time per pixel acquisition (s)

### 3.1.5 Detector Parameters

- **DetectorType:** Type of EBIC detector used
- **DetectorGain:** Detector gain setting

### 3.1.6 Image Acquisition Parameters

- **ImageWidth, ImageHeight:** Image dimensions in pixels
- **BitDepth:** Bit depth of the acquired image (typically 16-bit)
- **FrameCount:** Number of frames averaged

### 3.1.7 Sample and Stage Parameters

- **SampleBias:** Applied bias voltage during measurement (V)
- **StagePositionX, StagePositionY, StagePositionZ:** Sample stage coordinates

### 3.1.8 Instrument Identification

- **InstrumentModel, InstrumentSerial:** SEM instrument identification

### 3.1.9 Calibration Parameters

- **CurrentSensitivity:** Current measurement sensitivity
- **VoltageRange:** Voltage measurement range

These parameters are essential for accurate current conversion, spatial calibration, and reproducible EBIC measurements. The metadata extraction ensures proper scaling from pixel values to physical units and provides complete documentation of measurement conditions.

## 3.2 Current Map Conversion (CurrentMap Class)

The current conversion process involves two main steps: first converting pixel values to voltage measurements, then applying the instrument-specific calibration equation to obtain current values.

### 3.2.1 Pixel to Voltage Conversion

Raw 16-bit pixel values (0-65535) are converted to voltage using a linear scaling relationship:

$$V_{\text{measured}} = \left( \frac{\text{pixels}}{65535} \right) \times \text{scale} + \text{offset} \quad (4)$$

where scale = 1 mV and offset = -0.5 mV are instrument-specific calibration constants.

### 3.2.2 Voltage to Current Conversion

The fundamental equation used for converting measured voltage to current is derived from the SEM instrument's signal chain calibration:

$$I = \frac{\left( \frac{V_{\text{measured}} - O}{C} \right) \pm I_{\text{offset}}}{G} \times 10^9 \quad (5)$$

```

1 class CurrentMap:
2     def _compute_current(self, pixels):
3         # Extract calibration parameters from metadata
4         C = self.metadata['Contrast']           # Contrast setting
5         G = self.metadata['EffectiveAmpGain']    # Amplifier gain
6         O = self.metadata['OutputOffset']        # Output offset voltage
7         I_offset = self.metadata['InputOffset']  # Input offset
8         correction
9         inv = self.metadata['InverseEnabled']    # Signal inversion flag
10
11         # Step 1: Convert 16-bit pixels to voltage (mV)
12         scale = 1      # mV full-scale range
13         offset = -0.5  # mV zero offset
14         voltage = (pixels / 65535) * scale + offset
15
16         # Step 2: Apply instrument calibration equation
17         if inv:
18             # For inverted signal: I = [((V - O)/C) + I_offset] / G *
19             -1e9 nA
20             return (((voltage - O) / C) + I_offset) / G * -1e9
21         else:
22             # For normal signal: I = [((V - O)/C) - I_offset] / G * +1
23             e9 nA
24             return (((voltage - O) / C) - I_offset) / G * +1e9

```

Listing 2: Current conversion algorithm

The conversion algorithm accounts for all instrument-specific parameters:

- **Contrast (C)**: Amplifies or attenuates the signal range
- **Effective Amplifier Gain (G)**: Overall gain of the signal chain
- **Output Offset (O)**: DC offset voltage correction

- **Input Offset (I.offset):** Additional offset correction term
- **Inversion Mode:** Accounts for signal polarity inversion

The final result is expressed in nanoamperes (nA) by multiplying by  $10^9$ , providing physically meaningful current values for EBIC analysis. This calibrated current map forms the basis for all subsequent junction detection and diffusion length calculations.

### 3.2.3 SEMProcessor Class

The `SEMProcessor` class serves as the central coordinator that integrates all analysis components and manages the complete EBIC processing workflow:

```
1 class SEMProcessor:
2     def __init__(self, config):
3         self.config = config
4         self.metadata_parser = Metadata()
5         self.current_converter = CurrentMap()
6         self.junction_detector = JunctionAnalyzer()
7         self.diffusion_extractor = DiffusionLengthExtractor()
8         self.multi_image_analyzer = MultiImageAnalyzer()
```

Listing 3: Main processor class orchestrating the analysis pipeline

The `SEMProcessor` class provides:

- **Unified Interface:** Single entry point for all EBIC analysis operations
- **Workflow Management:** Coordinates the sequential execution of metadata extraction, current conversion, junction detection, and diffusion analysis
- **Error Handling:** Implements robust error handling and validation across all processing stages
- **Batch Processing Support:** Handles both single-image analysis and multi-image bias series
- **Result Aggregation:** Compiles results from all components into comprehensive output reports

This orchestrator class ensures that all analysis components work together seamlessly, providing a streamlined workflow from raw TIFF data to final EBIC parameter extraction.



### 3.3 Junction Detection Algorithms

#### 3.3.1 JunctionAnalyzer Class Overview

The `JunctionAnalyzer` class is a crucial component of the EBIC analysis framework. Its main role is to automatically and accurately detect the p-n junction line within a processed EBIC map. The class integrates filtering, edge detection, line fitting, coordinate transformation, and validation into a single robust workflow.

It contains methods for:

- **Preprocessing:** Filtering to reduce noise while keeping junction edges sharp.
- **Junction Detection:** Locating the junction using Canny edge detection and gradient-based analysis.
- **Post-processing:** Refining the noisy detection points into a smooth line using fitting techniques.
- **Coordinate Mapping:** Translating local ROI detections into global image coordinates.
- **Evaluation:** Quantitatively comparing the automatic detection with the manual reference.

#### 3.3.2 Preprocessing with Bilateral Filtering

The preprocessing stage applies a Bilateral Filter to the region of interest (ROI).

- **Purpose:** To reduce random noise while preserving the sharp junction edge.
- **Mechanism:** The bilateral filter considers both pixel intensity and spatial distance, smoothing flat areas without blurring steep gradients.

```
1 def _apply_preprocessing_filter(self, roi):
2     """
3     Applies a Bilateral Filter to smooth noise while preserving sharp
4     edges.
5     """
6     # Convert the ROI to 8-bit for OpenCV processing
7     normalized_roi = cv2.normalize(roi, None, 0, 255, cv2.NORM_MINMAX).
8     astype(np.uint8)
9
10    # Apply bilateral filter with optimized parameters
11    return cv2.bilateralFilter(normalized_roi, d=9, sigmaColor=75,
12                               sigmaSpace=9)
```

Listing 4: Bilateral filter preprocessing

#### 3.3.3 Canny Edge Detection with Otsu Thresholding

The core junction detection relies on column-wise edge analysis.

- **Adaptive Thresholding:** Otsu's method automatically determines thresholds for Canny.

- **Canny Edge Detection:** Extracts strong edge features while suppressing weak or false edges.
- **Gradient-Based Selection:** Ensures that the selected edge corresponds to the steepest EBIC transition.

```
1 def _detect_junction_canny(self, roi):
2     """
3     Detect junction using Canny edge detection with Otsu's adaptive
4     thresholding.
5     """
6     H, W = roi.shape
7     detected_points = np.zeros((W, 2), dtype=float)
8
9     # Convert ROI to 8-bit for OpenCV processing
10    roi_8bit = roi.astype(np.uint8)
11
12    # Otsu's method for automatic threshold determination
13    otsu_val, _ = cv2.threshold(roi_8bit, 0, 255, cv2.THRESH_BINARY +
14    cv2.THRESH_OTSU)
15
16    # Set Canny thresholds based on Otsu's result
17    high_threshold = float(otsu_val)
18    low_threshold = 0.5 * high_threshold
19
20    # Perform Canny edge detection
21    edges = cv2.Canny(roi_8bit, low_threshold, high_threshold)
22    .
23    .
24    .
25
26    return detected_points
```

Listing 5: Canny edge detection

### 3.3.4 Post-Processing for Junction Modeling

The detected points are refined into a continuous model.

**Linear Fitting** Applies linear regression to approximate the junction with a straight line.

```
1 def _fit_line_postprocessing(self, points):
2     """
3     Post-process detected points by fitting a linear model.
4     """
5     x_coords = points[:, 0]
6     y_coords = points[:, 1]
7
8     coeffs = np.polyfit(x_coords, y_coords, deg=1)
9     fitted_line = np.polyval(coeffs, x_coords)
10
11    return np.column_stack((x_coords, fitted_line))
```

Listing 6: Line fitting for junction approximation

**Spline Fitting** Uses cubic splines for more flexible modeling of curved junctions. **This method is not used for planar junctions. This is a good idea for curved junctions which we don't want to lose their structure while fitting**

```
1 def _fit_spline_postprocessing(self, points, smoothing=1.0):
2     """
3     Post-process detected points by fitting a spline model.
4     """
5     x_coords = points[:, 0]
6     y_coords = points[:, 1]
7
8     spline = UnivariateSpline(x_coords, y_coords, s=smoothing)
9     fitted_spline = spline(x_coords)
10
11     return np.column_stack((x_coords, fitted_spline))
```

Listing 7: Spline fitting for flexible junction modeling

### 3.3.5 Coordinate Transformation

The ROI-based detection is translated into global image coordinates.

```
1 def _map_detected_to_image_coords(self, detected_points, roi_coords):
2     """
3     Map detected points from ROI-local to global image coordinates.
4     """
5     x_offset, y_offset = roi_coords
6     mapped_points = detected_points.copy()
7     mapped_points[:, 0] += x_offset
8     mapped_points[:, 1] += y_offset
9     return mapped_points
```

Listing 8: Mapping detected junction to image coordinates

### 3.3.6 Validation Metrics

The automatic detection is validated against manual junctions using quantitative error metrics.

```
1 def _compare_with_manual(self, detected_points, manual_points):
2     """
3     Compare automatic and manual junction detection results.
4     """
5
6     return {
7         "mean_error": mean_error,
8         "max_error": max_error,
9         "rmse": rmse
10    }
```

Listing 9: Junction detection validation

## 3.4 Diffusion Length Extraction Algorithms

### 3.4.1 DiffusionLengthExtractor Class Overview

The `DiffusionLengthExtractor` class is a core component of the SEM-EBIC analysis software, designed to extract key physical parameters like diffusion length and depletion width from EBIC line profiles. The class's strength lies in its robust fitting algorithm, which is specifically engineered to handle noisy, real-world data and produce stable, reliable results through an innovative iterative fitting approach.

#### Class Purpose and Key Innovations

- **Iterative Stabilization:** Multi-pass fitting until parameter convergence
- **Head Search Algorithm:** Finding optimal starting points away from junction influence
- **Tail Truncation:** Systematic removal of noise-dominated regions
- **Convergence Detection:** Automated stability assessment for reliable results
- **Physical Parameter Extraction:** Diffusion length and depletion width calculation

### 3.4.2 Iterative Fitting Architecture

The core algorithm implements a sophisticated iterative process that addresses fundamental challenges in EBIC data analysis:

1. **Stage 1 - Head Search:** Progressive window shifting to find stable exponential decay regions
2. **Stage 2 - Tail Truncation:** Systematic removal of noisy tail regions
3. **Convergence Loop:** Repeated application until parameter stabilization
4. **Final Parameter Extraction:** Physical quantity calculation from stabilized fits

### 3.4.3 Stage 1: Head Search for Stable Starting Point

The `fit_profile_sides` method implements the critical first stage of finding robust starting points for exponential fitting.

```
1 def fit_profile_sides(self, x_vals, y_vals, intersection_idx,
2   max_shifts=20):
3     """
4     Stage 1: Find stable starting point through progressive window
5     shifting.
6     """
7
8     return stable_left_index, stable_right_index, left_results,
9     right_results
```

Listing 10: Head search with progressive window shifting

## Head Search Methodology

- **Progressive Shifting:** Moves fitting window incrementally away from junction peak
- **Exponential Model:** Fits standard minority carrier diffusion equation
- **Stability Monitoring:** Tracks diffusion length changes between shifts
- **Convergence Criterion:** Stops when changes fall below pixel-size tolerance
- **Physical Rationale:** Ensures fitting captures bulk material behavior, free from junction region influence

## Key Parameters

- `tolerance = self.pixel_size`: Convergence threshold based on physical resolution
- `max_shifts = 20`: Maximum iterations to prevent infinite loops
- `min_points = 10`: Minimum data points required for reliable fitting

### 3.4.4 Stage 2: Tail Truncation for Noise Mitigation

The `_truncate_tail_and_fit` method refines the fit by systematically removing noisy tail regions.

```
1 def _truncate_tail_and_fit(self, x_vals, y_vals, stable_start_idx, side
2     = 'left'):
3     """
4     Stage 2: Refine fit by progressively truncating noisy tail region.
5     """
6     except (RuntimeError, ValueError):
7         continue
8
9     # Select optimal result balancing fit quality and data utilization
10    if tail_results:
11        best_idx = np.argmax([r['r_squared'] * (0.3 + 0.7*r['
12        truncation_ratio'])
13        for r in tail_results])
14    return tail_results[best_idx]
15    else:
16        return None
```

Listing 11: Tail truncation for noise reduction

## Tail Truncation Strategy

- **Progressive Reduction:** Systematically shortens data from 100% to 50%
- **Quality Assessment:** Uses  $R^2$  values to evaluate fit goodness
- **Balanced Selection:** Considers both statistical quality and data utilization
- **Noise Elimination:** Removes low-SNR regions that distort exponential fits

### 3.4.5 Iterative Convergence Loop

The `fit_profile_sides_iterative` method orchestrates the complete iterative process.

```
1 def fit_profile_sides_iterative(self, x_vals, y_vals, intersection_idx,
2   max_iterations=5):
3     """
4     Main iterative fitting loop until parameter convergence.
5     """
6     prev_left_length = None
7     prev_right_length = None
8     convergence_tolerance = 1e-3 # Relative tolerance for convergence
9
10    return {
11        'left_diffusion_length': left_history[-1] if left_history else
12        None,
13        'right_diffusion_length': right_history[-1] if right_history
14        else None,
15        'depletion_width': depletion_width,
16        'iterations_used': iteration + 1,
17        'convergence_history': (left_history, right_history)
18    }
```

Listing 12: Iterative convergence algorithm

### Convergence Algorithm Features

- **Iterative Refinement:** Repeated application of head-tail process
- **Relative Tolerance:** Uses relative changes for convergence detection
- **History Tracking:** Maintains complete iteration history for analysis
- **Robust Handling:** Continues processing despite individual iteration failures
- **Convergence Detection:** Automatic stopping when parameters stabilize

### Convergence Criteria

- **Tolerance Setting:** `convergence_tolerance = 1e-3` (0.1% relative change)
- **Maximum Iterations:** `max_iterations = 5` prevents infinite loops
- **Bilateral Convergence:** Requires both left and right sides to stabilize
- **Physical Validation:** Ensures results represent true material properties

### 3.4.6 Physical Model and Parameter Extraction

**Exponential Decay Model** The core physical model describing minority carrier diffusion away from the junction:

$$I(x) = A \cdot \exp\left(-\frac{|x - x_j|}{L}\right) + B \quad (6)$$

Where:

- $A$ : Amplitude representing maximum EBIC current at junction
- $L$ : Diffusion length (characteristic decay distance)
- $x_j$ : Junction position
- $B$ : Background offset accounting for dark current

```
1 def _exponential_model(self, x, a, b, c):
2     """
3     Exponential decay model: y = a * exp(-b * x) + c
4     Corresponds to physical equation: I(x) = A * exp(-x/L) + B
5     where L = 1/b (diffusion length)
6     """
7     return a * np.exp(-b * x) + c
8
9 def _calculate_depletion_width(self, x_vals, left_index, right_index):
10    """
11    Calculate depletion width from stable fitting boundaries.
12    """
13    if left_index is None or right_index is None:
14        return None
15
16    left_boundary = x_vals[left_index]
17    right_boundary = x_vals[right_index]
18    depletion_width = abs(right_boundary - left_boundary) * self.
19    pixel_size_m
20
21    return depletion_width
```

Listing 13: Exponential model implementation

### 3.4.7 Parameter Initialization for Fitting

A critical part of the fitting process is providing good initial guesses for the fitting parameters. The `_fit_falling` helper function, which handles the actual curve fitting, intelligently sets these initial guesses (`p0`) based on the input data to ensure a robust and fast convergence. The `scipy.optimize.curve_fit` function then uses these guesses to find the best-fit parameters.

The exponential model being fitted is:

$$y = A \cdot e^{-\lambda(x-x_0)} + y_0$$

#### Initialization Strategy

- $y_0$  (**Offset**): Initialized as the last value of the profile, representing the background signal.
- $A$  (**Amplitude**): Estimated as the difference between the first data point (peak) and the last data point (baseline).
- $\lambda$  (**Decay Constant**): Set as the reciprocal of the total range of the x-values, providing a scale-aware starting guess.

- $x_0$  (**Shift**): Initialized as the first x-value in the data, corresponding to the junction peak.

This initialization scheme ensures fast and stable convergence even in the presence of noisy experimental EBIC profiles.

## 3.5 Interactive GUI for Visualization and Analysis

### 3.5.1 Interactive GUI (SEMViewer Class)

The GUI provides comprehensive visualization and interaction capabilities:

```
1 def _on_click(self, event): # Line drawing
2 def _on_scroll(self, event): # Zoom functionality
3 def _show_perpendicular_input(self, event): # Profile generation
4 def _fit_profiles(self, event): # Analysis triggering
```

Listing 14: GUI interaction handlers

Key features include:

- Multi-frame navigation between SEM and current maps
- Interactive line drawing for profile extraction
- Zoom and pan functionality with scroll and right-click drag
- Real-time perpendicular profile generation
- Junction detection visualization
- Diffusion length fitting and results display

## 3.6 Iteration over Multiple Images for Bias Analysis

### 3.6.1 Overview

The software implements a sophisticated pipeline for aligning multiple Scanning Electron Microscope (SEM) images and mapping key features such as p-n junctions and perpendicular profiles across different measurement conditions. This capability is essential for analyzing voltage-dependent effects and ensuring consistent region-of-interest tracking.

### 3.6.2 Global Image Alignment

**map\_all\_pixels: Multi-Dataset Alignment Orchestration** This high-level function manages the complete alignment process for multiple datasets:

```
1 def map_all_pixels(self, dataset_dict, window_size, max_disp, ref):
2     # Histogram matching for contrast normalization
3     matched_dict = {}
4     for key, (pix, cur) in dataset_dict.items():
5         matched_pix[0] = match_histograms(pix[0], image_ref)
6         matched_dict[key] = (matched_pix, cur)
7
8     # Compute mapping for each dataset
9     for key, (pix, cur) in matched_dict.items():
10        mapping = self.map_pixels(image_ref, pix[0],
```



```
11         window_size=window_size ,
12         search_radius=max_disp)
13     mapping_save[key] = mapping
```

Listing 15: Multi-dataset alignment orchestration

- **Histogram Matching:** Normalizes brightness and contrast variations across different measurements
- **Batch Processing:** Iterates through each dataset for comprehensive alignment
- **Result Storage:** Saves alignment mappings for subsequent analysis

**map\_pixels: Pixel-Level Correspondence** The core alignment algorithm implements robust pixel matching using Structural Similarity Index (SSIM):

```
1 def map_pixels(img1, img2, window_size=11, search_radius=15):
2     for i in range(height):
3         for j in range(width):
4             # Extract reference window
5             window1 = img1_padded[pi-pad:pi+pad+1, pj-pad:pj+pad+1]
6
7             # Search for best match in target image
8             for x in range(min_i, max_i):
9                 for y in range(min_j, max_j):
10                     window2 = img2_padded[x-pad:x+pad+1, y-pad:y+pad+1]
11                     score = ssim(window1, window2, data_range=
12
13                     data_range)
14
15                     if score > best_score:
16                         best_score = score
17                         best_pos = (x-pad, y-pad)
```

Listing 16: SSIM-based pixel matching

- **Structural Similarity Index (SSIM):** Advanced similarity metrics considering luminance, contrast, and structure
- **Brute-Force Search:** Comprehensive search within defined radius for optimal matches
- **Padding Handling:** Uses reflection padding to handle image boundaries

### 3.6.3 Line Alignment

**match\_line\_only: Optimized Feature Mapping** This function provides computational efficiency for line-based feature tracking:

```
1 def match_line_only(self, image1, image2, line_coords, step=5):
2     # Sample key points along the line
3     sampled_idx = np.arange(0, len(line_coords), step)
4
5     # Match only sampled pixels
6     for idx in sampled_idx:
7         x, y = line_coords[idx].astype(int)
8         # ... SSIM matching logic ...
```

```
9         mapping_sampled[idx] = best_pos
10
11     # Interpolate intermediate points
12     f_row = interp1d(idxs, rows, kind="linear")
13     f_col = interp1d(idxs, cols, kind="linear")
```

Listing 17: Optimized line matching with interpolation

- **Step-Wise Matching:** Reduces computation by matching key points only
- **Linear Interpolation:** Maintains line continuity between matched points
- **Performance Optimization:** Significant speed improvement for long lines

**match\_junction\_and\_perps: Structured Feature Mapping** Organizes the mapping of junction lines and perpendicular profiles:

```
1 def match_junction_and_perps(self, image1, image2, junction_line,
2   perpendiculars):
3     mapping_junction = self.match_line_only(image1, image2,
4     junction_line)
5
6     mapping_perps = []
7     for perp in perpendiculars:
8         m = self.match_line_only(image1, image2, perp)
9         mapping_perps.append(m)
```

Listing 18: Structured feature mapping

### 3.6.4 Final Data Extraction and Validation

**set\_line\_from\_mapping: Comprehensive Profile Extraction** The final workflow step generates physically meaningful data from aligned images:

```
1 def set_line_from_mapping(self, junction_line, mapping_junction, perps=
2   None, mapping_perps=None):
3     # Map junction line with fallback mechanisms
4     for (x, y) in junction_line.astype(int):
5         key = (int(y), int(x))
6         if key in mapping_junction:
7             r2, c2 = mapping_junction[key]
8             mapped_junc_list.append([c2, r2])
9         else:
10            # Neighborhood search fallback
11            for dy in (-1, 0, 1):
12                for dx in (-1, 0, 1):
13                    k = (int(y)+dy, int(x)+dx)
14                    if k in mapping_junction:
15                        r2, c2 = mapping_junction[k]
16                        mapped_junc_list.append([c2, r2])
17
18            # Interpolate missing points
19            valid = ~np.isnan(mapped_junction[:, 0])
20            if valid.sum() >= 2:
21                mapped_junction[:, 0] = np.interp(idxs, idxs[valid],
22                mapped_junction[valid, 0])
```

Listing 19: Profile extraction and validation

```
1 # Process each perpendicular profile
2 for i, (orig_perp, m) in enumerate(zip(perps, mapping_perps)):
3     # Map and interpolate points
4     # ...
5
6     # Calculate physical distances
7     dist_um = np.linspace(-((length_px-1)*self.pixel_size*1e6)/2,
8                           ((length_px-1)*self.pixel_size*1e6)/2,
9                           length_px)
10
11     # Sample SEM and current data
12     sem_vals = map_coordinates(sem_map, [ys, xs], order=1, mode="
nearest")
13     cur_vals = map_coordinates(current_map, [ys, xs], order=1, mode="
nearest")
```

Listing 20: Perpendicular profile processing

## Key Features of Data Extraction

- **Robust Mapping:** Fallback mechanisms handle imperfect matches through neighborhood searching
- **NaN Handling:** Intelligent interpolation maintains data continuity
- **Physical Unit Conversion:** Converts pixel distances to micrometers for quantitative analysis
- **Data Sampling:** Uses `map_coordinates` for precise SEM and EBIC data extraction
- **Intersection Detection:** Automatically finds junction intersection points

```
1 if show:
2     fig, ax = plt.subplots()
3     ax.imshow(self.pixel_maps[0], cmap="gray", origin="upper")
4     if self.detected_junction_line is not None:
5         ax.plot(self.detected_junction_line[:, 0],
6                 self.detected_junction_line[:, 1],
7                 color="green", linewidth=2, label="Mapped junction")
8     # Plot perpendiculars and intersection points
```

Listing 21: Visual validation

This multi-image alignment capability enables precise comparative analysis of semiconductor properties under different bias conditions, providing crucial insights into voltage-dependent behavior of diffusion lengths and depletion regions.

## **4 Analysis Workflow**

### **4.1 Step 1: Data Loading and Conversion**

1. User selects TIFF files through GUI
2. Software extracts XML metadata and parses instrument parameters
3. Converts pixel values to current maps using calibration data
4. Saves processed data as CSV files for each frame

### **4.2 Step 2: Junction Detection**

1. User draws approximate line along suspected junction
2. Software extracts rectangular ROI perpendicular to the line
3. Applies Canny edge detection with spline smoothing
4. Returns refined junction position with accuracy metrics

### **4.3 Step 3: Profile Analysis**

1. Generates perpendicular profiles across the junction
2. Extracts SEM contrast and current values along each profile
3. Marks intersection points with detected junction
4. Applies low-pass filtering to reduce high-frequency noise

### **4.4 Step 4: Exponential Fitting**

1. Fits falling exponentials to both sides of each profile
2. Uses shifted fitting windows to achieve stability
3. Calculates inverse diffusion lengths ( $1/\lambda$ ) for each fit
4. Computes average depletion widths and diffusion lengths
5. Generates quality metrics ( $R^2$  values) for each fit

## 5 Software Workflow Overview

This section describes the comprehensive workflow of the SEM-EBIC analysis tool. Currently, two distinct versions are available: a GUI-focused implementation emphasizing user interaction and visualization, and a characterization-oriented version prioritizing analytical capabilities. The long-term objective involves merging these implementations into a unified codebase while extending the overall functionality.

### 5.1 GUI-Based Software Workflow

#### 5.1.1 TIFF File Processing

The analysis workflow commences with the processing of TIFF image files. The initial interface presents a file browser enabling users to select directories containing SEM data in TIFF format.




Figure 1: Directory selection interface for TIFF file processing

Upon directory selection, the system displays all available TIFF files within the chosen folder. Users may select individual files or multiple files for batch processing. Subsequently, the software prompts for an output directory where processed results will be stored.




Figure 2: TIFF file selection interface showing available datasets

Following processing completion, the interface presents a list of available images for visualization and analysis. Users can select specific datasets for detailed examination.




Figure 3: Processed file selection interface for data visualization

### 5.1.2 Main Analysis Interface

The primary analysis window provides comprehensive access to all analytical tools and functionalities. This centralized interface facilitates seamless navigation through the complete EBIC analysis pipeline.




Figure 4: Main analysis interface with complete tool accessibility

The interface includes navigation controls for multi-frame datasets, allowing sequential examination of different image frames and toggling between various data representations.




Figure 5: Frame navigation and display mode toggle controls

### 5.1.3 Image Visualization Tools

The software incorporates advanced visualization capabilities, including overlay functionality for comparative analysis of different data channels or processing results.




Figure 6: Data overlay functionality for multi-channel visualization

The zoom functionality employs a rectangular selection mechanism. Users activate the zoom tool, delineate the region of interest, and subsequently deactivate the tool to return to full-image viewing.




Figure 7: Interactive zoom feature with rectangular region selection

#### 5.1.4 Junction Analysis Workflow

The junction detection process initiates with manual line drawing approximating the p-n junction location. The system subsequently prompts for window width specification to define the analysis region, optimizing junction detection while mitigating interference from extraneous features.






Figure 8: Manual junction approximation with adjustable analysis window

The perpendicular profile analysis module generates cross-sectional profiles across the detected junction. Users specify the quantity and length of perpendicular lines, with the resulting profiles displaying both SEM intensity and EBIC current data. Intersection points between profiles and the junction are explicitly marked at the zero-position reference.




Figure 9: Perpendicular profile analysis showing SEM and EBIC data intersections




Figure 10: Detected junction and manual junction with perpendicular line




Figure 11: Depletion region characterization and width quantification

A visualization of the data used for fitting is overlaid on the original profile, providing an intuitive understanding of the fitting process. The visualization shows different truncations on tail which is shortening the tail. The left side is inverted to suit falling exponential model.




Figure 12: Right tail truncation




Figure 13: Left tail truncation

## 5.2 Parameter Extraction and Results

### 5.2.1 Iterative Stabilization Algorithm for Diffusion Length Extraction

The core innovation in the diffusion length extraction methodology is the iterative stabilization algorithm that ensures robust and physically meaningful parameter estimation. The process follows a sophisticated multi-stage approach:

#### 5.2.2 Algorithm Workflow

The stabilization process operates through the following sequential steps:

1. **Initial Shift Phase:** Systematic shifting of the fitting window starting point both left and right away from the junction position
2. **Exponential Fitting:** Application of nonlinear least-squares fitting to the exponential decay model  $I(x) = A \cdot e^{-\lambda(x-x_0)} + y_0$  for each shifted window position
3. **Convergence Monitoring:** Continuous evaluation of the extracted decay constant  $\lambda$  (and thus diffusion length  $L = 1/\lambda$ ) until stability is achieved
4. **Tail Truncation Phase:** Progressive removal of data points from the tail region where signal-to-noise ratio deteriorates
5. **Stability Verification:** Repeated fitting with truncated data to confirm parameter stability
6. **Iterative Refinement:** Cyclic repetition of the entire process until both sides of the junction demonstrate convergence

### 5.2.3 Convergence Criteria Enhancement

As noted, the convergence condition can be improved by requiring bilateral agreement between the left and right side fittings. The enhanced convergence criterion becomes:

$$\text{Converged} = \left( |L_{\text{left}}^{(k)} - L_{\text{left}}^{(k-1)}| < \epsilon \right) \wedge \left( |L_{\text{right}}^{(k)} - L_{\text{right}}^{(k-1)}| < \epsilon \right) \wedge \left( |L_{\text{left}}^{(k)} - L_{\text{right}}^{(k)}| < \delta \right) \quad (7)$$

Where  $\delta$  represents an additional tolerance for inter-side consistency, ensuring that the extracted parameters are physically plausible and consistent across both sides of the junction.

Average Lengths Summary

Metric	Average Value
Depletion Width ( $\mu\text{m}$ )	0.126
Diffusion Length Left ( $\mu\text{m}$ )	0.288
Diffusion Length Right ( $\mu\text{m}$ )	0.048
Diffusion Length (All) ( $\mu\text{m}$ )	0.216

Pixel size = 0.006  $\mu\text{m}$

Figure 14: Comprehensive fitting results with quantitative parameters

The depletion region analysis provides specific characterization of the junction properties, displaying the spatially resolved depletion width measurements.




Figure 15: Fitted profiles at 0 original starting point (0 shift)

### 5.2.4 Enhancement of Physical Model: Bi-exponential Diffusion

The current model is a simple exponential is a simple exponential as discussed before, while it can be improved in the way that is discussed below.

Minority carrier concentration decay with distance is modeled using a bi-exponential function:

$$C(x) = A_1 e^{-x/L_1} + A_2 e^{-x/L_2} \quad (8)$$

where  $L_1$ ,  $L_2$  are diffusion lengths of two transport components, and  $A_1$ ,  $A_2$  are their relative amplitudes.

The two terms reflect distinct mechanisms:

- $L_1, A_1$ : bulk diffusion with longer characteristic length
- $L_2, A_2$ : surface or defect-related recombination with shorter length

A signal-weighted average provides a practical single parameter:

$$L_{\text{eff}} = A_1 L_1 + A_2 L_2 \quad (9)$$

### 5.2.5 Significance

Compared to single-exponential models, the bi-exponential approach offers:

- better representation of multiple recombination pathways
- improved fit to complex decay profiles
- insight into both bulk and surface effects

Thus, while  $L_{\text{eff}}$  is a convenient summary, the full model captures the richer physics of carrier transport.

## 5.3 Pixel matching

You can go ahead with drawing junction and defining perpendicular profiles, the red line is the manual line and the green line is the detected junction.

keep in mind that this version is using Spline post processing which means that the detected junction is not a line. 3.3.4

After completing the same workflow as done before, you can close every plot and also the first image, then processing of the second image starts. You should then see the last line in this image, meaning that the second image is started to be processed.

```
Profile 1: converged after 2 iterations.  
Figure closing, saving images and data...  
Saved overlay image to tiff_test_output\sample1\saved_on_close\sample1_frame1_overlay.png  
Warning: pixel_map.csv not found in output_root\sample2\frame_2
```

Figure 16: Second image processing

Then the matched junction and perpendicular line show up, as it's obvious the length and orientation has changed a little. Change of orientation is ok, but the length should not change. The change in the length is caused by line fitting that may include pixels that has moved alot.




Figure 17: matched lines

Then you can fit profiles as you did before, the results, somehow, match our predictions.




Figure 18: Depletion region defined after pixel matching

### 5.3.1 Profile Matching Precision

Since the detected profile is not perfectly perpendicular to the junction, the fitting is not sufficiently precise. This limitation can be mitigated by improving the algorithm used for matching the profiles across images. Also simple line fitting is not the final sufficient destiny, the procedure can be improved. The method is noted here 3.6.3




Figure 19: Bad profiles for fitting

## **6 Future Development Plans**

### **6.1 Enhanced Edge Detection Based on Current Map**

Development of advanced edge detection algorithms specifically optimized for EBIC current maps, improving junction detection accuracy and robustness.

### **6.2 Software Version Integration**

Integration of the currently separate GUI-based and characterization-focused versions into a unified codebase with comprehensive functionality.

### **6.3 Convergence Criteria Enhancement**

As discussed in Section 5.2.3

#### **6.3.1 Enhancement of Physical Model: Bi-exponential Diffusion**

As discussed in Section 5.2.4

### **6.4 Advanced Pixel Matching Algorithms**

Implementation of more sophisticated pixel matching techniques for improved multi-image alignment and feature tracking across different measurement conditions.

### **6.5 Pixel Matching and Profile Matching Precision**

As discussed in Section 5.3.1, the profile alignment can be improved ...